

# EIS Analysis Toolkit

**Version:** v0.9.4 (2026-01-05)

Modular toolkit for electrochemical impedance spectroscopy (EIS) analysis with Distribution of Relaxation Times (DRT) support.

## Key features:

- **Reproducibility** - Objective, data-driven parameter selection
- **Modularity** - Usable as CLI and Python library
- **Advanced optimization** - Multi-start, Differential Evolution

## Supported data formats:

- **Gamry DTA** - native format (automatic metadata parsing, ZCURVE block)
- **CSV** - three columns with header: frequency, Z\_real, Z\_imag
  - Delimiter: comma, semicolon, or tab (auto-detection)
  - Decimal format: US (dot) and European (comma for semicolon-delimited)
  - Comments: lines starting with # are ignored
  - Examples: [example/example\\_eis\\_data.csv](#)

**Changelog:** [CHANGELOG.md](#)

---

## Quick start

### Installation

```
git clone https://github.com/chmelat/eis_analysis
cd eis_analysis
pip install -e .
```

After installation, the `eis` command is available system-wide.

For alternative installation methods, see [Installation options](#).

### Basic usage

```
# Basic analysis (KK validation + DRT)
eis data.DTA

# With circuit fitting
eis data.DTA --circuit "R(100) - (R(5000) | C(1e-6))"

# Save plots to files
eis data.DTA --save results --format pdf
```

Run `eis --help` for all options.

---

## Common workflows

### Data quality check

```
# KK validation + DRT analysis (default)
eis data.DTA
```

Output: Kramers-Kronig residuals, mu metric, DRT spectrum with detected peaks.

### Circuit fitting

```
# Fit equivalent circuit
eis data.DTA --circuit "R(100) - (R(5000) | C(1e-6))"

# Use Differential Evolution for global optimization (default)
eis data.DTA --circuit "R(100) - (R(5000) | Q(1e-6, 0.9))"

# Use multi-start for local optimization
eis data.DTA --circuit "..." --optimizer multistart --multistart 20
```

### Automatic circuit suggestion

```
# Automatic Voigt chain fitting
eis data.DTA --voigt-chain

# With automatic element count optimization
eis data.DTA --voigt-chain --voigt-auto-M
```

### Oxide layer analysis

```
# Thickness calculation from capacitance
eis data.DTA --circuit "R(100) - (R(5000) | C(1e-6))" --analyze-oxide

# With custom parameters
eis data.DTA --circuit "..." --analyze-oxide --epsilon-r 22 --area 0.5
```

### Batch processing

```
# Process multiple files without interactive display
for f in *.DTA; do
    eis "$f" --save "${f%.DTA}" --no-show
done
```

---

## Main features

### Kramers-Kronig validation

Data quality verification using Lin-KK test. Validates causality, linearity, and stability of measured data.

*# Default (included in standard analysis)*

```
eis data.DTA
```

*# Skip KK validation*

```
eis data.DTA --no-kk
```

*# Custom mu threshold*

```
eis data.DTA --mu-threshold 0.80
```

**Detailed documentation:** [doc/LinKK\\_analysis.md](#)

## DRT analysis

Distribution of Relaxation Times - model-free method for impedance data analysis. Automatic regularization parameter selection using GCV (Generalized Cross-Validation).

*# Default DRT with auto-lambda*

```
eis data.DTA
```

*# Manual lambda selection*

```
eis data.DTA --lambda 1e-3
```

*# GMM peak detection (more robust)*

```
eis data.DTA --peak-method gmm
```

*# Robust R<sub>inf</sub> estimation for inductive data*

```
eis data.DTA --ri-fit
```

*# Peak classification into physical processes*

```
eis data.DTA --peak-method gmm --classify-terms
```

**Detailed documentation:** [doc/GCV\\_IMPLEMENTATION.md](#), [doc/GMM\\_PEAK\\_DETECTION.md](#)

## Circuit fitting

Elegant operator overloading syntax for circuit definition.

### Supported elements:

Element	Description	Example
R(value)	Resistor	R(100)
C(value)	Capacitor	C(1e-6)
L(value)	Inductor	L(1e-6)
Q(Q, n)	Constant Phase Element (CPE)	Q(1e-4, 0.8)
W(sigma)	Warburg (semi-infinite)	W(50)
Wo(R <sub>W</sub> , tau)	Warburg (bounded)	Wo(100, 1.0)
K(R, tau)	Voigt with tau parametrization	K(1000, 1e-4)

Values in parentheses serve as initial guesses for the nonlinear fitting algorithm. Values in quotes (e.g., `R("100")`) are treated as fixed constants and will not be fitted.

### Operators:

Operator	Meaning	Example
-	Series connection	<code>R(100) - C(1e-6)</code>
	Parallel connection	<code>R(1000)   C(1e-6)</code>

**Operator precedence:** - has HIGHER precedence than | (Python rules). Always use parentheses around parallel combinations: `(R|C)`.

*# Voigt element*

```
eis data.DTA --circuit "R(100) - (R(5000) | C(1e-6))"
```

*# Randles circuit with CPE*

```
eis data.DTA --circuit "R(10) - (R(100) | Q(1e-4, 0.8))"
```

*# With fixed parameter*

```
eis data.DTA --circuit 'R("0.86") - (R(2.4e9) | Q(1e-10, 0.823))'
```

**Detailed documentation:** [doc/CIRCUIT\\_PARSER.md](#), [doc/K\\_ELEMENT\\_GUIDE.md](#)

### Voigt chain (automatic circuit)

Automatic Voigt chain estimation using linear regression for initial guess, then nonlinear refinement.

*# Basic Voigt chain*

```
eis data.DTA --voigt-chain
```

*# Automatic element count optimization*

```
eis data.DTA --voigt-chain --voigt-auto-M
```

*# Custom density (elements per decade)*

```
eis data.DTA --voigt-chain --voigt-n-per-decade 5
```

**Detailed documentation:** [doc/VOIGT\\_CHAIN\\_MATH.md](#)

### Advanced optimization

**Differential Evolution (default)** Global optimization for finding global minimum.

*# Default DE*

```
eis data.DTA --circuit "R(100) - (R(5000) | C(1e-6))"
```

*# Custom parameters*

```
eis data.DTA --circuit "..." --de-strategy 2 --de-popsize 20 --de-maxiter 500
```

**DE strategies:** 1=randtobest1bin (default), 2=best1bin, 3=rand1bin

**Detailed documentation:** [doc/DIFFERENTIAL\\_EVOLUTION.md](#)

**Multi-start optimization** Multiple starts from different initial points.

*# Multi-start with 20 restarts*

```
eis data.DTA --circuit "... " --optimizer multistart --multistart 20
```

*# With larger perturbation*

```
eis data.DTA --circuit "... " --optimizer multistart --multistart-scale 3.0
```

**Detailed documentation:** [doc/MULTISTART\\_OPTIMIZATION.md](#)

## Oxide layer analysis

Oxide layer thickness calculation from capacitance.

*# Automatic analysis*

```
eis data.DTA --circuit "... " --analyze-oxide
```

*# Custom permittivity and area*

```
eis data.DTA --circuit "... " --analyze-oxide --epsilon-r 9 --area 0.5
```

Common permittivities: ZrO<sub>2</sub> ~ 22, Al<sub>2</sub>O<sub>3</sub> ~ 9, TiO<sub>2</sub> ~ 80, SiO<sub>2</sub> ~ 3.9

**Detailed documentation:** [doc/OXIDE\\_ANALYSIS\\_GUIDE.md](#)

---

## CLI Reference

### Input data and frequency range

- input - Input file (.DTA or .csv). Without argument, synthetic data is used for testing. Built-in test circuit: **Rs - (R0||Q0) - (R1||Q1)** with 1% Gaussian noise:
  - Rs = 10  $\Omega$  (series resistance)
  - R0 = 100 k $\Omega$ , Q0 = (1e-6 S·s<sup>n</sup>, n=0.6)
  - R1 = 800 k $\Omega$ , Q1 = (3e-5 S·s<sup>n</sup>, n=0.43)
- --f-min - Minimum frequency [Hz]. Data below this value will be cut off. Useful for removing noise at low frequencies.
- --f-max - Maximum frequency [Hz]. Data above this value will be cut off. Useful for removing artifacts at high frequencies.

### Circuit fitting

- --circuit, -c - Equivalent circuit for fitting. Syntax: - = series, | = parallel. Example: "R(100) - (R(5000) | C(1e-6))". Supported elements: R, C, L, Q, W, Wo, K.
- --weighting (default: sqrt) - Weighting type for fitting: uniform (all points equal), sqrt (square root of 1/|Z|), proportional (1/|Z|), square (1/|Z|^2).
- --no-fit - Skip circuit fitting.

## Optimizer selection

- `--optimizer` (default: `de`) - Optimizer type: `de` (Differential Evolution - global), `multistart` (multiple local fits), or `single` (one local fit).

## Differential Evolution (global optimization)

- `--de-strategy` (default: `1`) - DE strategy: `1=randtobest1bin` (balanced, default), `2=best1bin` (fast convergence), `3=rand1bin` (more exploration).
- `--de-popsiz` (default: `15`) - Population size as multiple of parameter count. Higher = better exploration but slower.
- `--de-maxiter` (default: `1000`) - Maximum number of generations. Increase if optimization doesn't converge.
- `--de-tol` (default: `0.01`) - Convergence tolerance (relative fitness change).
- `--de-workers` (default: `1`) - Number of parallel workers. `-1` = all CPU cores.

## Multi-start optimization

- `--multistart N` - Number of restarts for multi-start optimization. Each restart starts from perturbed initial values. Default: `16` when `--optimizer multistart` is used.
- `--multistart-scale` (default: `2.0`) - Perturbation size in sigma units (standard deviation from covariance matrix). Higher = larger parameter space exploration.

## DRT analysis

- `--lambda, -l` (default: `auto GCV`) - Regularization parameter for DRT. Without this parameter, automatic selection using GCV (Generalized Cross-Validation) and L-curve method is used. Higher values = smoother DRT, lower = more detail but also noise.
- `--n-tau, -n` (default: `100`) - Number of points on the tau time constant axis. Higher values give finer DRT resolution but increase computational cost.
- `--normalize-rpol` - Normalize  $\gamma(\tau)$  by polarization resistance so that  $\int \gamma(\tau) d\tau = 1$ . Useful for comparing samples with different  $R_{pol}$ .
- `--peak-method` (default: `scipy`) - Peak detection method in DRT: `scipy` (fast, `scipy.signal.find_peaks`) or `gmm` (robust, Gaussian Mixture Model - requires `scikit-learn`).
- `--gmm-bic-threshold` (default: `10.0`) - BIC threshold for GMM peak detection. Lower values detect more peaks (`2-5`: sensitive, `10-20`: conservative). Only used with `--peak-method gmm`.
- `--ri-fit` - Robust  $R_{inf}$  estimation using  $R+L+K$  model fit on high-frequency data. Suitable for data with inductive loop.
- `--classify-terms` - Classification of DRT peaks into physical processes (requires `--peak-method gmm`).
- `--no-drt` - Skip DRT analysis. Useful if you only want circuit fitting.

## Kramers-Kronig validation

- `--no-kk` - Skip Kramers-Kronig validation. KK test verifies causality, linearity, and stability of data.
- `--mu-threshold` (default: 0.85) - Threshold value of mu metric for Lin-KK test. Values below threshold indicate problematic data.
- `--ocv` - Display OCV (Open Circuit Voltage) curve if available in data.

## Voigt chain (automatic circuit)

- `--voigt-chain` - Use automatic Voigt chain fitting. Linear regression for  $R_i$  and  $\tau_i$  estimation, then nonlinear refinement.
- `--voigt-n-per-decade` (default: 3) - Number of time constants per decade for Voigt chain. Higher = finer coverage but more parameters.
- `--voigt-extend-decades` (default: 0.0) - Extend tau range by N decades beyond data limits. Useful if you expect processes outside measured range.
- `--voigt-prune-threshold` (default: 0.01) - Threshold for removing small  $R_i$  (as fraction of  $R_{pol}$ ). Elements with  $R_i < \text{threshold} * R_{pol}$  are removed.
- `--voigt-allow-negative` - Allow negative  $R_i$  values (Lin-KK style). Otherwise negative elements are removed.
- `--voigt-no-inductance` - Do not include series inductance L in model.
- `--voigt-fit-type` (default: complex) - Fit type: complex (default, real+imag), real (real part only), imag (imaginary part only).
- `--voigt-auto-M` - Automatically optimize number of M elements using mu metric.
- `--voigt-mu-threshold` (default: 0.85) - Mu threshold value for `--voigt-auto-M`.
- `--voigt-max-M` (default: 50) - Maximum number of M elements for `--voigt-auto-M`.
- `--no-voigt-info` - Do not display detailed Voigt chain fitting info.

## Oxide layer analysis

- `--analyze-oxide` - Perform oxide layer analysis - thickness calculation from capacitance.
- `--epsilon-r` (default: 22.0) - Relative permittivity of oxide. Default 22 for ZrO<sub>2</sub>. Other oxides: Al<sub>2</sub>O<sub>3</sub> ~ 9, TiO<sub>2</sub> ~ 80, SiO<sub>2</sub> ~ 3.9.
- `--area` (default: 1.0) - Electrode area [cm<sup>2</sup>]. Required for correct thickness calculation.

## Jacobian

- `--numeric-jacobian` - Use numeric Jacobian instead of analytic. Analytic Jacobian is faster and more accurate but not available for all elements. Use this option for custom/non-standard elements.

## Output and visualization

- `--save, -s` - Save plots with this prefix. Example: `--save results` creates `results_nyquist.png`, etc.

- `--format, -f` (default: png) - Format of saved plots: png (raster), pdf/svg/eps (vector for publications).
- `--no-show` - Do not display plots interactively. Useful for batch processing with `--save`.
- `-v, --verbose` - Show debug messages on stderr (prefix [DEBUG]).
- `-q, --quiet` - Quiet mode - hide INFO messages, show only warnings and errors.

### Logging levels:

- INFO (default): stdout, no prefix
  - WARNING: stdout, prefix !
  - ERROR: stderr, prefix !!
  - DEBUG (`-v`): stderr, prefix [DEBUG]
- 

## Installation options

### Option 1: pip install (recommended)

For an isolated environment, we recommend using a virtual environment (venv). This prevents dependency conflicts with other projects and does not affect your system Python installation. If you prefer not to use venv, skip to the direct installation below.

#### With virtual environment (recommended):

Linux/macOS:

```
# Create virtual environment
python3 -m venv eis_env

# Activate environment
source eis_env/bin/activate

# Install (now in isolated environment)
git clone https://github.com/chmelat/eis_analysis
cd eis_analysis
pip install -e .      # Editable install (for development)
# or
pip install .         # Standard install
```

Windows:

```
# Create virtual environment
python -m venv eis_env

# Activate environment (may require execution policy, see below)
eis_env\Scripts\Activate.ps1

# Install
git clone https://github.com/chmelat/eis_analysis
cd eis_analysis
pip install -e .
```



For Windows, you may need to set the execution policy:

```
Set-ExecutionPolicy -ExecutionPolicy RemoteSigned -Scope CurrentUser
```

### Without virtual environment (direct installation):

```
git clone https://github.com/chmelat/eis_analysis
cd eis_analysis
pip install -e .      # Editable install (for development)
# or
pip install .         # Standard install
```

After installation, the eis command is available (in activated environment if using venv):

```
eis --help
eis data.DTA
```

### Option 2: Run without installation

If you prefer not to install the package, you can run the script directly:

```
pip install numpy scipy matplotlib      # Install dependencies
python3 eis.py --help                  # Run script directly
python3 eis data.DTA
```

### Option 3: System packages (Debian/Ubuntu)

```
sudo apt install python3-numpy python3-scipy python3-matplotlib
pip install -e .      # Then install the package
```

### Optional dependencies

```
# GMM peak detection
pip install scikit-learn

# Development tools (ruff, mypy, pytest)
pip install -e ".[dev]"
```

### Windows

```
git clone https://github.com/chmelat/eis_analysis
cd eis_analysis
pip install -e .
eis data.DTA
```

---

## Python API

Complete Python API: [doc/PYTHON\\_API.md](doc/PYTHON_API.md)

---

## Documentation

Document	Description
<a href="#">doc/PYTHON_API.md</a>	Complete Python API reference
<a href="#">doc/CIRCUIT_PARSER.md</a>	Circuit parser syntax
<a href="#">doc/K_ELEMENT_GUIDE.md</a>	K element guide
<a href="#">doc/LinKK_analysis.md</a>	Kramers-Kronig validation
<a href="#">doc/GCV_IMPLEMENTATION.md</a>	GCV technical documentation
<a href="#">doc/GMM_PEAK_DETECTION.md</a>	GMM peak detection
<a href="#">doc/DRT_METHOD_ANALYSIS.md</a>	DRT method analysis
<a href="#">doc/VOIGT_CHAIN_MATH.md</a>	Voigt chain mathematics
<a href="#">doc/DIFFERENTIAL_EVOLUTION.md</a>	Differential Evolution
<a href="#">doc/MULTISTART_OPTIMIZATION.md</a>	Multi-start optimization
<a href="#">doc/NONLINEAR_FIT_ANALYSIS.md</a>	Nonlinear optimization overview
<a href="#">doc/OXIDE_ANALYSIS_GUIDE.md</a>	Oxide layer analysis
<a href="#">doc/CODE_ANALYSIS_REPORT.md</a>	Code structure analysis

---

## References

- Orazem, M.E., Tribollet, B.: *Electrochemical Impedance Spectroscopy* (2008)
- Schonleber, M. et al.: "A Method for Improving the Robustness of linear Kramers-Kronig Validity Tests", *Electrochimica Acta* 131 (2014)
- Wahba, G.: "A comparison of GCV and GML", *Annals of Statistics* 13 (1985)
- Saccoccio, M. et al.: "Optimal regularization in DRT", *Electrochimica Acta* 147 (2014)

---

This code was developed with the assistance of [Claude Code](#).

**License:** MIT License